



World Anti-Doping Program

## GUIDELINES

# **Human GROWTH HORMONE (hGH) BIOMARKERS TEST**

*for Doping Control Analyses*

Version 2.0

April 2016

## TABLE OF CONTENTS

<b>1. Objective</b>	<b>3</b>
<b>2. Scope</b>	<b>3</b>
<b>3. Introduction to the Method</b>	<b>3</b>
3.1 Principle of the Method	4
<b>4. Assay Requirements</b>	<b>5</b>
4.1 Assay Pre-analytical Procedure	6
4.2 Assay Analytical Procedure	9
4.2.1 Analytical Testing Strategy	9
<b>5. Reporting and Interpretation of Results</b>	<b>11</b>
5.1 Interpretation of Test Results	11
5.1.1 Presumptive Adverse Analytical Finding (PAAF)	13
5.1.2 Adverse Analytical Finding (AAF)	13
5.1.3 Atypical Finding (ATF)	13
5.2 Reporting of Test Results	14
<b>6. Assay Measurement Uncertainty</b>	<b>15</b>
6.1 Combined Standard Uncertainty ( $u_c$ )	15
6.2 Maximum levels of $u_c$	16
6.3 Expanded Uncertainty ( $U_{95\%}$ )	16
6.4 Verification of Measurement Uncertainty	16
<b>7. Definitions</b>	<b>17</b>
7.1 Code Defined Terms	17
7.2 ISL Defined Terms	18
7.3 International Standard for Testing and Investigations (ISTI) Defined Terms	19
<b>8. Bibliography</b>	<b>20</b>
Extended Bibliography describing the hGH Biomarkers Test	22

## **1. Objective**

This guideline has been developed to ensure a harmonized approach in the application of the GH-2000 Biomarkers Test for the detection of doping with human Growth Hormone (hGH) in sport. The guideline provides direction on the *Sample* pre-analytical preparation procedure, the performance of the test and the interpretation of the test results.

## **2. Scope**

This guideline follows the rules established in the *World Anti-Doping Agency's (WADA) International Standards for Laboratories (ISL)* [1] and relevant Technical Documents regarding the testing of blood *Samples*. These requirements are still fully applicable and shall be respected. This guideline contains additional recommendations to facilitate the implementation of the *Testing* procedures particular to hGH detection by the Biomarkers Test.

## **3. Introduction to the Method**

The hGH Biomarkers Test involves the measurement of two hGH-sensitive *Markers*, namely insulin-like growth factor-I (IGF-I) and N-terminal pro-peptide of type III collagen (P-III-NP), which are present in serum. The Bibliography at the end of these guidelines lists the main publications produced during the development and validation of the method. These measurements are combined in sex-specific discriminant function formulae which improve the sensitivity and specificity of the test based on a score (the GH-2000 score) [2] to detect hGH misuse compared with single-*Marker* analysis. The hGH Biomarkers Test may also have utility in detecting GH secretagogues and IGF-I abuse in sport [3, 4].

A series of placebo-controlled recombinant (r)hGH administration studies performed in Europe (lead centers in the UK and Germany) and Australia has shown that both IGF-I and P-III-NP rise substantially following rhGH administration in a dose-dependent manner [2, 5-11]. These *Markers* have been evaluated for several confounding factors that might influence the scores of the discriminant functions, including age, sex [2], ethnicity [12], exercise [8, 9], diurnal and day-to-day variation, intra-individual variation [13], bony and soft tissue injury [14], sporting discipline, and body habitus (physique) [15-17].

Except sex and age, no other factor has been shown to affect the hGH discriminant function scores substantially.

The GH-2000 discriminant function formulae are sex-specific, based on the natural logarithm of IGF-I and P-III-NP serum concentrations (required to normalize the data distribution) and include an adjustment for age to reflect the age-related decline in hGH and *Marker* concentrations [2].

### 3.1 Principle of the Method

The hGH Biomarkers Test is based on the measurement of IGF-I and P-III-NP by immunoassays or mass spectrometry (MS)-based approaches [18].

In order to perform the test, an assay pairing formed by an IGF-I and a P-III-NP assay is utilized for the Initial Testing Procedure, whereas two different IGF-I/P-III-NP assay pairings shall be used for the Confirmation Procedures (see Table 2 below). One IGF-I/P-III-NP assay pairing may be the same as that used in the Initial Testing Procedure. It is recommended that the Liquid Chromatography (LC)-tandem MS (LC-MS/MS) or LC-High Resolution MS (LC-HRMS) assay for IGF-I be applied as part of the Confirmation Procedure whenever possible. The results of each assay pairing are then used to calculate the GH-2000 score.

The assays currently used are:

#### IGF-I assays

- 1) Immunotech A15729 IGF-I IRMA assay (Immunotech SAS, Marseille, France)

The Immunotech assay is a two-site, solid-phase, immunoradiometric assay (IRMA) using two monoclonal antibodies prepared against two different antigenic sites of the IGF-I molecule. The first is coated on a solid phase and the second is radiolabelled with  $^{125}\text{I}$ . IGF-I is separated from IGFBPs by acidification and excess IGF-II is added to prevent further interference with the assay from IGFBPs. The Immunotech assay is calibrated using the WHO IGF-I IRP standard 87/518.

- 2) IDS-iSYS IGF-I assay (Immunodiagnosics Systems Limited, Boldon, UK).

The iSYS IGF-I assay is an automated sandwich, chemiluminescent immunoassay (CLIA). Samples are incubated with an acidic solution to dissociate IGF-I from the IGFBPs. A portion of this, along with a neutralization buffer containing excess IGF-II to prevent re-aggregation with IGFBPs, a biotinylated anti-IGF-I monoclonal antibody directed against the N-terminal, and an acridinium labeled anti-IGF-I monoclonal antibody are incubated. Streptavidin labeled magnetic particles are then added and, following an additional incubation step, the magnetic particles are captured using a magnet. After a washing step and addition of trigger reagents, the light emitted by the acridinium label is directly proportional to the concentration of IGF-I in the original sample [19]. The iSYS IGF-I assay is calibrated using the new WHO recombinant IGF-I IRP standard 02/254.

- 3) LC-MS/MS or LC-HMRS IGF-I assay [18].

This is a bottom-up approach based on the quantification of peptides derived from trypsin digestion of IGF-I. Serum samples are incubated with an acidic solution in the presence of excess IGF-II and  $^{15}\text{N}$ -labeled IGF-I as internal standard. Proteins are precipitated with acetonitrile. Following reduction and alkylation of the dried supernatant, the solution is enzymatically hydrolyzed with trypsin. Two peptides corresponding to amino acids 1–21 (T1) and 22–36 (T2) of IGF-I are separated by LC and measured by MS/MS or HRMS.

## P-III-NP assays

### 1) Orion UniQ™ P-III-NP RIA (Orion Diagnostica, Espoo, Finland)

The Orion UniQ™ P-III-NP RIA is a competitive radioimmunoassay based on the formation of a complex between solid-phase anti-P-III-NP polyclonal rabbit antibodies and P-III-NP in the serum samples in competition with <sup>125</sup>I-labelled P-III-NP. A sample volume of 100 µL is used.

### 2) Siemens ADVIA Centaur P-III-NP assay [(Siemens Healthcare Laboratory Diagnostics, Camberley, UK)] [20]

The Siemens ADVIA Centaur P-III-NP assay is an automated, two-site sandwich, chemiluminescent immunoassay. The assay uses two monoclonal mouse antibodies: the first antibody is an acridinium ester-labeled anti-P-III-NP antibody. The second antibody is a biotin-labeled anti-P-III-NP antibody. The solid phase contains streptavidin-coated paramagnetic particles and during the reaction, the light emitted by the acridinium label is directly proportional to the concentration of P-III-NP in the sample. The Siemens P-III-NP assay is calibrated by the manufacturer using a standard derived from bovine P-III-NP.

## **4. Assay Requirements**

Prior to the implementation of the Biomarkers Test in routine *Doping Control* analysis, the Laboratory shall fulfill the following requisites:

- Validate the assays' performance on-site, including the determination of the assays' Limit of Quantification (*LOQ*), Repeatability ( $s_r$ ), Intermediate Precision ( $s_w$ ) and bias;
- The acceptance values for parameters of assay performance, applicable to the separate determinations of IGF-I and P-III-NP concentrations, are specified in Table 1 below;
- In addition, the Laboratory shall determine the assay Measurement Uncertainty (*MU*) from Laboratory validation data. The combined standard uncertainty ( $u_c$ ) shall be not higher than a maximum value of  $u_{c\_Max} = 0.50$  for either assay pairing, expressed as Standard Deviations (*SD*) and applied to the GH-2000 scores at values close to the corresponding Decision Limits (*DLs*), as described in section 7 below;
- Demonstrate readiness for assay implementation through test validation data and/or successful participation in at least one *WADA*-approved educational External Quality Assessment Scheme (*EQAS*) round or inter-Laboratory collaborative study. In cases of identified deficiencies, proper corrective action(s) shall be documented and implemented;
- Obtain ISO/IEC-17025 accreditation from a relevant accreditation body for the inclusion of the hGH Biomarkers Test in the Laboratory scope of accreditation.

**Table 1:** Acceptance Criteria for parameters of assay performance.

Validation parameter	Immunoassays	LC-MS/MS or LC-HRMS <sup>a</sup>
<b><math>s_r</math></b> (within-assay Relative Standard Deviation, <i>RSD</i> %)	≤ 10%	≤ 10%
<b><math>s_w</math></b> (between-assay <i>RSD</i> %)	≤ 20%	≤ 15%
<b>LOQ<sup>b</sup></b> IGF-I P-III-NP	≤ 50 ng/mL ≤ 1 ng/mL	≤ 50 ng/mL N/A

<sup>a</sup> when applied to the mean of the measured concentrations of T1 and T2.

<sup>b</sup> LOQ is defined as the lowest concentration meeting the criteria for  $s_r$  and  $s_w$ .

#### 4.1 Assay Pre-analytical Procedure

Upon reception of the "A" and "B" *Samples* in the Laboratory, the following steps should be followed:

- Check that the blood *Samples* have been collected in tubes containing an inert polymeric serum separator gel and a clotting activation factor (BD Vacutainer® SST™-II Plus tubes, EU ref 367955; BD Vacutainer® SST™-II Plus *Advance* tubes, EU ref 367954) in accordance with the *WADA* Guidelines for Blood *Sample* Collection [21]. Such blood *Samples* should have been kept in a refrigerated state (not frozen) following collection and during transportation to the Laboratory<sup>1</sup>;
- Alternatively, *Samples* may be received in the Laboratory as frozen or refrigerated serum *Samples*, following the clotting and centrifugation of the blood and separation of the serum fraction at the site of *Sample* collection;
- Any *Samples* delivered to the Laboratory as plasma shall not be accepted for the purposes of hGH analysis with the current assays. In line with this, the Sample Collection Authorities are provided with Guidelines for collection of blood *Samples* for hGH analysis, which specify that the matrix of analysis is serum [21]. The Laboratory shall notify and seek advice from the Testing Authority regarding rejection and Analytical Testing of *Samples* for which irregularities are noted (as per ISL 6.2.2.4). In cases of *Sample* collection in the incorrect matrix (to be identified at the results management level), the results of such analysis of the *Sample* shall be disregarded;

<sup>1</sup> Previous studies have demonstrated that IGF-I and P-III-NP concentrations remain stable if the sample remains refrigerated for up to 5 days [22].

- Check the status of the *Sample(s)* (e.g. evidence of hemolysis) and the integrity of the collection tubes (e.g. evidence of breakage of the separating gel). The Laboratory shall note any unusual condition of the *Sample*, record such condition(s) and include it in the Test Report to the Testing Authority;
- For *Samples* received as whole blood in SST™-II tubes or SST™-II Plus Advance tubes:

#### "A" *Sample*

- Centrifuge the "A" *Sample* for 10-15 min at 1300-1500 g as soon as possible after reception;
- The whole separated serum fraction from the "A" *Sample* should be transferred into another tube or aliquoted into new vials, which shall be properly labelled to ensure Laboratory Internal Chain of Custody documentation. One Aliquot should be used for the Initial Testing Procedure. The remaining "A" *Sample Aliquot(s)* not used for the Initial Testing Procedure must be stored frozen<sup>2</sup> until the "A" Confirmation Procedure, if needed;
- For the Initial Testing Procedure, "A" *Sample Aliquots* may be analyzed immediately after aliquoting or stored at approximately 4°C for a maximum of 24h before analysis (within a maximum of 5 days from *Sample* collection). Alternatively, the "A" *Sample Aliquots* must be frozen<sup>2</sup> until analysis.

#### "B" *Sample*

- Centrifuge the "B" *Sample* for 10-15 minutes at 1300-1500g as soon as possible after reception. The whole of the "B" *Sample* separated serum fraction should be kept in the SST™-II or SST™-II Plus Advance *Sample*

---

<sup>2</sup> For storage of Aliquots frozen, well-closing vials should be used (for optimal storage cryovials with an "O-ring" are recommended) and the following conditions are recommended:

- For short-term storage (up to three months) at approximately -20°C;
- For long-term periods (more than three months) freeze at approximately -20°C and transfer to approximately -70 to -80°C.

Thawing of the *Sample(s)* for analysis shall not be done under hot water or any other similar process that would raise the temperature of the *Sample* above room temperature. Thawing overnight at approximately 4°C is recommended.

collection tube and step-frozen (refrigeration prior to freezing) according to the tube manufacturer's instructions<sup>3</sup> until analysis, if needed;

- Once the "B" *Sample* is thawed and opened (according to ISL 6.2.4.2.2), an Aliquot of the "B" *Sample* shall be used for the "B" Confirmation Procedure. The remaining "B" *Sample* serum should be transferred into a new tube/vial and shall be sealed in front of the *Athlete* or the *Athlete's* representative or a Laboratory-appointed independent witness using a tamper-proof evident method and frozen<sup>2</sup> until further analysis, if needed.
- For *Samples* received as separated serum *Samples*:
  - a) *Samples* received as frozen separated serum fractions:
    - These *Samples* should remain frozen<sup>2</sup> until analysis;
    - Once thawed, an Aliquot of *Sample* "A" shall be taken to be used for the Initial Testing Procedure. This Aliquot of *Sample* "A" may be stored at approximately 4°C if the Initial Testing Procedure is scheduled to take place within 24h of thawing. The remaining "A" *Sample* serum fraction may be kept in the *Sample* collection tube or aliquoted into new vials, which shall be properly labelled to ensure Laboratory Internal Chain of Custody documentation, and stored frozen<sup>2</sup> until the "A" Confirmation Procedure, if needed;
    - Once the "B" *Sample* is thawed and opened (according to ISL 6.2.4.2.2), an Aliquot of the "B" *Sample* shall be used for the "B" Confirmation Procedure. The remaining "B" *Sample* serum shall be kept in the *Sample* collection tube and shall be sealed in front of the *Athlete* or the *Athlete's* representative or a Laboratory-appointed independent witness using a tamper-proof evident method and frozen<sup>2</sup> until further analysis, if needed.
  - b) *Samples* received as refrigerated separated serum fractions:
    - Take an Aliquot of the "A" *Sample* as soon as possible upon reception. For the Initial Testing Procedure, "A" *Sample Aliquots* may be analyzed immediately after aliquoting or stored at approximately 4°C for a maximum of 24h before analysis (within a maximum of 5 days from *Sample* collection). Alternatively, "A" *Sample Aliquots* must be frozen<sup>2</sup> until analysis;
    - The remainder of the "A" *Sample* not used for the Initial Testing Procedure may be kept in the *Sample* collection tube or aliquoted into new vials,

---

<sup>3</sup> Place the tube into a dedicated isolating box before transferring into a -20°C freezer. In order to maintain the integrity of the separation gel, allow the freezing to proceed for at least 2 hours before moving or transferring the frozen tubes. Moving the tubes before the separating gel is frozen and stable may lead to contamination of serum by cellular material.



which shall be properly labelled to ensure Laboratory Internal Chain of Custody documentation, and stored frozen<sup>2</sup> until the "A" Confirmation Procedure, if needed;

- For "B" *Samples*, freeze<sup>2</sup> the *Samples* as soon as possible upon reception and thaw before analysis. Once the "B" *Sample* is thawed and opened (according to ISL 6.2.4.2.2), an Aliquot of the "B" *Sample* shall be used for the "B" Confirmation Procedure. The remaining "B" *Sample* serum shall be kept in the *Sample* collection tube and shall be re-sealed in front of the *Athlete* or the *Athlete's* representative or a Laboratory-appointed independent witness using a tamper-proof evident method and stored frozen<sup>2</sup> until further analysis, if needed.

## 4.2 Assay Analytical Procedure

- For the performance of the assay(s) analytical procedure, refer to the test procedure described in the Instructional Insert provided with the test assays and the Laboratory SOP;
- In cases of contradiction between the Instructional Insert provided with the assays and the Laboratory SOP, or between the Instructional Insert and these Guidelines, the latter document shall prevail in each case.

### 4.2.1 Analytical Testing Strategy

- One assay pairing (e.g. Immunosoft IGF-I + Orion P-III-NP) should be used for the Initial Testing Procedure (Table 2);
- In the case of an initial Presumptive Adverse Analytical Finding (PAAF), two different assay pairings shall be used for the Confirmation Procedure of the "A" *Sample* (Table 2) using three new Aliquots of the original "A" *Sample*<sup>4</sup>. One of the assay pairings may be the same as the one used for the Initial Testing Procedure;
- For the "B" Confirmation Procedure, both assay pairings used during the confirmation of the "A" *Sample* shall be applied on three Aliquots taken from the original "B" *Sample*<sup>5</sup>. The Laboratory shall follow the requirements of the ISL 6.2.4.2.2.1 for the performance of the "B" *Sample* confirmation analysis;

---

<sup>4</sup> Laboratories that do not have the analytical capacity to perform the Confirmation Procedure with an additional assay pairing shall have, upon consultation with the responsible Testing Authority, the *Sample* shipped to and analyzed by another Laboratory that has such analytical capacity.

- Either the LC-MS/MS or LC-HRMS IGF-I assay may be applied as the unique test for IGF-I quantification (*i.e.* either assay may be used for the Initial Testing Procedure and also combined with two different P-III-NP assays for the Confirmation Procedure(s));
- For both "A" and "B" Confirmation Procedures, three Sample Aliquots shall be measured, except in cases of limited Sample volume, in which case a lower maximum number of replicates may be used;
- In accordance with the ISL provisions 6.2.4.2.1.4 and 6.2.4.2.2.8, the Laboratory shall have a policy to define those circumstances where the Confirmation Procedure of an "A" or "B" Sample should be repeated (for example, values of within-assay  $s_r > 10\%$ );
- It is recommended that the Laboratories implement well-characterized and stable internal quality control (QC) sample(s), which are not subject to assay lot variations, for the performance of the tests under different assay conditions (different lots of assay, different analysts, *etc.*). Following preparation/reception by the Laboratory, all QC material should be aliquoted and stored frozen (preferably at  $-80^\circ\text{C}$  for long-term storage) until use.

These QC samples<sup>5</sup> should be:

- **QC<sub>low</sub>**: Serum obtained from healthy individual(s), which is shown to have a value of  $\leq 200$  ng/mL IGF-I and  $< 5$  ng/mL P-III-NP;
- **QC<sub>high</sub>**: Serum obtained from hGH administration studies or another appropriate source that has been shown to contain concentrations of  $\geq 500$  ng/mL IGF-I and  $\geq 10$  ng/mL P-III-NP.
- Assay Repeatability ( $s_r$ ) and Intermediate Precision ( $s_w$ ) will be assessed by analyzing each QC sample in triplicates on 5-6 separate occasions;
- With every new batch of reagents (new lot number), the following evaluation steps should be implemented before accepting the new batch:
  - Each of the QC samples shall be determined at least three times whenever a new batch of reagents is obtained. The number of replicates per determination shall be as stipulated by the assay manufacturers. The QCs may be measured in a single assay or over a range of assays. If, for any QC, the difference between the mean concentration for the new batch and that for the preceding batch is more than 20%, investigation of the new batch will be required;

---

<sup>5</sup> Four QC samples may also be used, as long as they contain IGF-I and P-III-NP at the necessary concentrations (*e.g.* QC<sub>IGF-I<sub>low</sub></sub>, QC<sub>IGF-I<sub>high</sub></sub>, QC<sub>P-III-NP<sub>low</sub></sub> and QC<sub>P-III-NP<sub>high</sub></sub>).

- In order to detect small but systematic changes with time, it is recommended that the performance of a new batch of reagents be controlled, for example, through a cumulative sum (CUSUM) chart/table, which is built for each QC based on the difference between the mean(s) for the new batch and the initial value(s). When using the CUSUM, results should be assessed using customary procedures as detailed at

<http://itl.nist.gov/div898/handbook/pmc/section3/pmc323.htm>;

## **5. Reporting and Interpretation of Results**

### **5.1 Interpretation of Test Results**

For determination of compliance of the analytical result, the Laboratory shall compare the *Sample's* GH-2000 score (rounded to two decimal places) with the corresponding gender-specific DLs established for the assay pairings [23].

- The DL values are given in Table 2 below <sup>6</sup>;
- The MU of the assays has already been considered and incorporated in the reference population-based statistical estimation of the DL<sup>7</sup> [24, 25]. Therefore, for declaration of an *AAF* or an *ATF* the assay MU shall not be added.

---

<sup>6</sup> The DL values specified above have been derived from the analysis of *Samples* from *Athletes* treated under *Doping Control* conditions of *Sample* collection, transportation, storage and analysis [23]. The established DL values define a combined test specificity (between the two assay pairings used for the Confirmation Procedure) of at least 99.99%. These DL values are conservative values and will be periodically refined as more data are accumulated from normative studies and *Doping Control* tests performed by WADA-accredited laboratories.

<sup>7</sup> According to WADA's Technical Document on Decision Limits for the Confirmatory Quantification of Threshold Substances (TDDL) [24], the decision rule applicable to assays for which the Threshold value(s) have been established based on reference population statistics already incorporates a guard band that reflects the uncertainty of the measurements provided by the assay(s). Therefore, the zone of analytical values considered compliant (negative) or not (*AAF*) with this decision rule would be defined by the Threshold value itself, which constitutes the DL.

The GH-2000 score for the *Sample* is calculated applying the following discriminant function formulae:

GH-2000 score for males:

$$-6.586 + 2.905 \cdot \ln(P\text{-III-NP}) + 2.100 \cdot \ln(IGF\text{-I}) - 101.737 / \text{age}$$

GH-2000 score for females:

$$-8.459 + 2.454 \cdot \ln(P\text{-III-NP}) + 2.195 \cdot \ln(IGF\text{-I}) - 73.666 / \text{age}$$

where  $\ln(P\text{-III-NP})$  and  $\ln(IGF\text{-I})$  are the natural logarithms (ln) of the mean concentration values (expressed in ng/mL) obtained from the measured replicates of the *Sample Aliquot* and age is rounded down to the nearest year <sup>8</sup>.

**Table 2.** Possible assay pairings for the Initial Testing Procedure and Confirmation Procedure(s) and applicable sex-specific Decision Limits.

<b>Sex</b>	<b>Assay Pair</b> (IGF-I + P-III-NP)	<b>DL1</b>
<b>Males</b>	LC-MS/MS or LC-HRMS + Orion	9.70
	LC-MS/MS or LC-HRMS + Siemens Advia Centaur	11.34
	IDS-Sys + Orion	9.00
	IDS-Sys + Siemens Advia Centaur	10.61
	ImmunoTech + Orion	9.98
	ImmunoTech + Siemens Advia Centaur	11.53
<b>Females</b>	LC-MS/MS or LC-HRMS + Orion	8.56
	LC-MS/MS or LC-HRMS + Siemens Advia Centaur	10.13
	IDS-Sys + Orion	7.79
	IDS-Sys + Siemens Advia Centaur	9.35
	ImmunoTech + Orion	8.62
	ImmunoTech + Siemens Advia Centaur	10.10

<sup>8</sup> For calculation of the GH-2000 scores, the natural logarithms (ln) of the mean concentrations (ng/mL) of IGF-I and P-III-NP shall be expressed to 3 decimal places. However, for compliance decisions (comparison to the assay pairing- and gender-specific DLs), the resulting GH-2000 score shall be rounded to two decimal places.

#### 5.1.1 Presumptive Adverse Analytical Finding (PAAF)

- The Initial Testing Procedure shall produce a PAAF for *Sample "A"* if the corresponding GH-2000 score (rounded to two decimal places) exceeds the sex-specific DL (Table 2) applicable for the assay pairing used for the screening procedure;
- When the LC-MS/MS or LC-HRMS method is used for IGF-I quantification during the Initial Testing Procedure, the test result shall be considered a PAAF if the GH-2000 score, calculated on the basis of the IGF-I concentration determined from the quantification of the T1 or the T2 diagnostic peptide (Table 3), exceeds the sex-specific DL applicable for the assay pairing used (Table 2).

#### 5.1.2 Adverse Analytical Finding (AAF)

- The Confirmation Procedure shall produce an AAF if the *Sample's* GH-2000 scores (rounded to two decimal places) exceed the sex-specific DLs (Table 2) established for the two assay pairings applied for the Confirmation Procedure;
- When the LC-MS/MS or LC-HRMS method is used for IGF-I quantification during the Confirmation Procedure, the test result shall be considered an AAF if:
  - the GH-2000 scores calculated on the basis of the **average** IGF-I concentration determined from the quantification of T1 and T2 exceed the sex-specific DLs established in Table 2 for the two assay pairings applied, and the T1- and T2-derived IGF-I concentrations do not differ by more than 20% (Table 3).

#### 5.1.3 Atypical Finding (ATF)

- The Confirmation Procedure shall produce an ATF if the GH-2000 scores (rounded to two decimal places) exceed the DL (Table 2) for only one of the two assay pairings employed for the Confirmation Procedure;
- When the LC-MS/MS or LC-HRMS method is used for IGF-I quantification during the Confirmation Procedure, the test result shall also be considered an ATF if:
  - the GH-2000 scores calculated on the basis of the **average** IGF-I concentration determined from the quantification of T1 and T2 exceed the sex-specific DLs established in Table 2, BUT
  - the IGF-I concentrations determined from the quantification of T1 and T2 differ by more than 20% (Table 3);
  - In such cases, the Laboratory shall repeat the LC-MS/MS or LC-HRMS analysis to verify the IGF-I T1, T2 concentration difference before reporting the finding.

**Table 3.** Examples of interpretation of tests findings when applying LC-MS/MS or LC-HRMS for IGF-I quantification.

Procedure	$\frac{ T_1 - T_2 }{\text{MEAN}(T_1; T_2)}$	GH-2000 score			Interpretation/ Reporting
		IGF-I (T1)	IGF-I (T2)	Mean IGF-I (T1, T2)	
<u>Initial Testing Procedure</u>	N/A	N/A	> DL	N/A	<b>PAAF</b>
		> DL	N/A	N/A	<b>PAAF</b>
<u>Confirmation Procedure</u>	≤ 0.2	> DL	> DL	> DL	<b>AAF</b>
		> DL	< DL	> DL < DL	<b>AAF</b> <i>Negative</i>
		< DL	> DL	> DL < DL	<b>AAF</b> <i>Negative</i>
		< DL	< DL	< DL	<i>Negative</i>
	> 0.2	> DL	> DL	> DL	<b>ATF</b>
		> DL	< DL	> DL < DL	<b>ATF</b> <i>Negative</i>
		< DL	> DL	> DL < DL	<b>ATF</b> <i>Negative</i>
		< DL	< DL	< DL	<i>Negative</i>

## 5.2 Reporting of Test Results

- When reporting an *AAF* or an *ATF*, the Laboratory Test Report shall include the mean GH-2000 scores from triplicate determinations (obtained during the Confirmatory Procedure) expressed to two decimal places, the values of the applicable DL as well as the combined standard uncertainty of the assay ( $u_c$ , expressed as *SD*) at values close to the DL as determined by the Laboratory;
- In addition, the Laboratory Documentation Package shall include the mean concentration values of IGF-I and P-III-NP from triplicate determinations (obtained during the Confirmatory Procedure, expressed to the nearest integer for IGF-I and two decimal places for P-III-NP) and the expanded MU equivalent to the 95% coverage interval ( $U_{95\%}$ ,  $k = 2$ ) for the value of the GH-2000 score for the *Sample*.

**Test Report Example** (e.g. for a *Sample* from a male *Athlete*):

The analysis of the *Sample* with the hGH Biomarkers Test has produced the following GH-2000 scores: 10.90 for assay pair '1' [IDS IGF-I + Centaur P-III-NP] and 9.90 for assay pair '2' [LC-MS/MS IGF-I + Orion P-III-NP], which are greater than the corresponding male-specific DLs of 10.61 and 9.70, respectively. The combined standard uncertainty ( $u_c$ ) estimated by the Laboratory at levels close to the DL is 0.40 for assay pair '1' and 0.35 for assay pair '2'. This constitutes an *Adverse Analytical Finding* for hGH.

## 6. Assay Measurement Uncertainty

### 6.1 Combined Standard Uncertainty ( $u_c$ )

- Laboratories shall generally refer to the TDDL [24] for estimation of assay MU;
- The Laboratories shall determine each assay's  $u_c$  based on their assay validation data;

The  $u_c$  is a dynamic parameter that can be reduced with increasing improvement in the performance of the assays. The establishment of a confident value of  $u_c$  would be based on multiple measurements done throughout a long period of time, when certain sources of uncertainty (such as environmental changes, instrument performance, different analysts, etc.) would be accounted for;

- ISO/IEC 17025 recommends that  $u_c$  be estimated using an approach consistent with the principles described in the ISO/IEC Guide to the Expression of Uncertainty in Measurement (GUM) [26];
- For application to the hGH marker method, the following approach for calculation of the  $u_c$  budget is recommended:

The value of  $u_c$ , applicable to the GH-2000 scores close to the DLs, will result from the contributing  $u_c$  of the component assays (applicable to the natural logarithms (ln) of the values of the measured concentrations) using the law of propagation of uncertainty, according to formulae (1)<sup>9</sup>:

$$(1) \text{ For males: } u_c (\text{score}) = \sqrt{8.44 * u_c^2 [\ln (\text{P-III-P})] + 4.41 * u_c^2 [\ln (\text{IGF-I})]}$$

$$\text{For females: } u_c (\text{score}) = \sqrt{6.02 * u_c^2 [\ln (\text{P-III-P})] + 4.82 * u_c^2 [\ln (\text{IGF-I})]}$$

---

<sup>9</sup> In formula (1) and (2), the  $u_c$  (score) and the contributing  $u_c$  associated with the values of the natural logarithms of the measured concentrations should be expressed as standard deviations (*SD*).

- The  $u_c$  associated with the values of the natural logarithms (ln) of the concentrations determined with the IGF-I and P-III-NP assays, shall be estimated from the Intermediate Precision ( $s_w$ ) and the bias of the ln determinations according to formula (2)<sup>9</sup>;

$$(2) \quad u_c = \sqrt{s_w^2 + u_{bias}^2}$$

- For calculation of  $u_c$ , either a single QC sample, containing IGF-I and P-III-NP in appropriate concentrations (e.g. QC<sub>high</sub>) or two separate QC samples containing IGF-I at ~500-800 ng/mL (e.g. QC<sub>IGFI-high</sub>) and P-III-NP at ~10-20 ng/mL (e.g. QC<sub>PIIIINP-high</sub>), should be used<sup>10</sup>. These QCs should be aliquoted and stored frozen (preferably at -80°C for long term storage) until use;
- QC sample(s) and four different ½ dilutions should be measured in triplicates over 5-6 days by at least 2 different analysts. This would ensure that the  $s_w$  is calculated over the physiological range of concentrations of hGH *Markers* that may be found in samples producing GH-2000 scores close to the DLs;
- The bias will be established by comparison of the Laboratory's long-term means of the ln of concentration values obtained e.g. for the QC<sub>low</sub> and QC<sub>high</sub> samples with the expected values determined through a *WADA* educational EQAS round or inter-Laboratory collaborative study. The bias contribution to  $u_c$  is expressed as  $RMS_{bias}$ .

## 6.2 Maximum levels of $u_c$

In accordance with the TDDL [24], Laboratories shall have values of  $u_c$ , applicable to values close to the DL for each assay pairing, not higher than the maximum values of  $u_{c\ Max}$ .

## 6.3 Expanded Uncertainty ( $U_{95\%}$ )

For determination of the expanded uncertainty  $U_{95\%}$  a coverage factor  $k=2$  can be applied if  $u_c$  has a 95 % confidence level.

$$(3) \quad U_{95\%} = k * u_c, \text{ where } k=2$$

## 6.4 Verification of Measurement Uncertainty

Laboratories shall refer to the TDDL [24] for ongoing verification of the assay MU estimates.

---

<sup>10</sup> Since the GH-2000 scores depend on the age of the donor, in order to produce relevant values of the GH-2000 scores (close to the DLs), the age of the donors should ideally be between 20 – 40 years old.



## **7. Definitions**

### **7.1 Code Defined Terms**

*Adverse Analytical Finding:* A report from a WADA-accredited laboratory or other WADA - approved laboratory that, consistent with the International Standard for Laboratories and related Technical Documents, identifies in a *Sample* the presence of a *Prohibited Substance* or its *Metabolites* or *Markers* (including elevated quantities of endogenous substances) or evidence of the *Use* of a *Prohibited Method*.

*Athlete:* Any *Person* who competes in sport at the international level (as defined by each International Federation) or the national level (as defined by each *National Anti-Doping Organization*). An *Anti-Doping Organization* has discretion to apply anti-doping rules to an *Athlete* who is neither an *International-Level Athlete* nor a *National-Level Athlete*, and thus to bring them within the definition of "Athlete." In relation to *Athletes* who are neither *International-Level* nor *National-Level Athletes*, an *Anti-Doping Organization* may elect to: conduct limited *Testing* or no *Testing* at all; analyze *Samples* for less than the full menu of *Prohibited Substances*; require limited or no whereabouts information; or not require advance *TUEs*. However, if an Article 2.1, 2.3 or 2.5 anti-doping rule violation is committed by any *Athlete* over whom an *Anti-Doping Organization* has authority who competes below the international or national level, then the *Consequences* set forth in the *Code* (except Article 14.3.2) must be applied. For purposes of Article 2.8 and Article 2.9 and for purposes of anti-doping information and education, any *Person* who participates in sport under the authority of any *Signatory*, government, or other sports organization accepting the *Code* is an *Athlete*.

*Atypical Finding:* A report from a WADA-accredited laboratory or other WADA -approved laboratory which requires further investigation as provided by the International Standard for Laboratories or related Technical Documents prior to the determination of an *Adverse Analytical Finding*.

*Code:* The World Anti-Doping Code.

*Doping Control:* All steps and processes from test distribution planning through to ultimate disposition of any appeal including all steps and processes in between such as provision of whereabouts information, *Sample* collection and handling, laboratory analysis, *TUEs*, results management and hearings.

*International Standard:* A standard adopted by WADA in support of the *Code*. Compliance with an *International Standard* (as opposed to another alternative standard, practice or procedure) shall be sufficient to conclude that the procedures addressed by the *International Standard* were performed properly. *International Standards* shall include any Technical Documents issued pursuant to the *International Standard*.

*Marker:* A compound, group of compounds or biological variable(s) that indicates the *Use* of a *Prohibited Substance* or *Prohibited Method*.

*Sample* or *Specimen:* Any biological material collected for the purposes of *Doping Control*.

*Testing:* The parts of the *Doping Control* process involving test distribution planning, *Sample* collection, *Sample* handling, and *Sample* transport to the laboratory.

*WADA:* The World Anti-Doping Agency.

## 7.2 ISL Defined Terms

Aliquot: A portion of the *Sample* of biological fluid or tissue (e.g. urine, blood) obtained from the *Athlete* used in the analytical process.

Analytical Testing: The parts of the *Doping Control* process involving *Sample* handling, analysis and reporting following receipt in the Laboratory.

Confirmation Procedure: An analytical test procedure whose purpose is to identify the presence or to measure the concentration/ratio of one or more specific *Prohibited Substances, Metabolite(s)* of a *Prohibited Substance*, or *Marker(s)* of the *Use of a Prohibited Substance* or *Method* in a *Sample*.

[*Comment: A Confirmation Procedure for a Threshold Substance shall also indicate a concentration/ratio of the *Prohibited Substance* greater than the applicable Decision Limit (as noted in the TD DL).*]

Decision Limit: a concentration, accounting for the maximum permitted combined uncertainty, above which an *Adverse Analytical Finding* shall be reported.

Initial Testing Procedure: An analytical test procedure whose purpose is to identify those *Samples* which may contain a *Prohibited Substance, Metabolite(s)* of a *Prohibited Substance*, or *Marker(s)* of the *Use of a Prohibited Substance* or *Prohibited Method* or the quantity of a *Prohibited Substance, Metabolite(s)* of a *Prohibited Substance*, or *Marker(s)* of the *Use of a Prohibited Substance* or *Prohibited Method*.

Intermediate Precision: Variation in results observed when one or more factors, such as time, equipment, or operator are varied within a Laboratory.

*International Standard* for Laboratories (ISL): The *International Standard* applicable to Laboratories as set forth herein.

Laboratory Internal Chain of Custody: Documentation of the sequence of *Persons* in custody of the *Sample* and any Aliquot of the *Sample* taken for Analytical Testing.

[*Comment: Laboratory Internal Chain of Custody is generally documented by a written record of the date, location, action taken, and the individual performing an action with a *Sample* or Aliquot.*]

Laboratory(ies): (A) WADA-accredited laboratory(ies) applying test methods and processes to provide evidentiary data for the detection of *Prohibited Substances, Methods* or *Markers* on the *Prohibited List* and, if applicable, quantification of a Threshold Substance in *Samples* of urine and other biological matrices in the context of anti-doping activities.

Laboratory Documentation Packages: The material produced by the Laboratory to support an analytical result such as an *Adverse Analytical Finding* as set forth in the WADA Technical Document for Laboratory Documentation Packages.

Measurement Uncertainty (MU): Parameter associated with a measurement result that characterizes the dispersion of quantity values attributed to a measurand.

[*Comment: Knowledge of the MU increases the confidence in the validity of a measurement result*].

Presumptive Adverse Analytical Finding: The status of a *Sample* test result for which there is a suspicious result in the Initial Testing Procedure, but for which a confirmation test has not yet been performed.

Repeatability, sr: Variability observed within a Laboratory, over a short time, using a single operator, item of equipment, etc.

Threshold Substance: An exogenous or endogenous *Prohibited Substance, Metabolite or Marker* of a *Prohibited Substance* which is analyzed quantitatively and for which an analytical result (concentration, ratio or score) in excess of a pre-determined Decision Limit constitutes an *Adverse Analytical Finding*. Threshold Substances are identified as such in the Technical Document on Decision Limits (TD DL).

### **7.3 International Standard for Testing and Investigations (ISTI) Defined Terms**

Sample Collection Authority: The organization that is responsible for the collection of *Samples* in compliance with the requirements of the *International Standard for Testing and Investigations*, whether (1) the Testing Authority itself; or (2) another organization (for example, a third party contractor) to whom the Testing Authority has delegated or sub-contracted such responsibility (provided that the Testing Authority always remains ultimately responsible under the *Code* for compliance with the requirements of the *International Standard for Testing and Investigations* relating to collection of *Samples*).

Testing Authority: The organization that has authorized a particular *Sample* collection, whether (1) an *Anti-Doping Organization* (for example, the International Olympic Committee or other *Major Event Organization, WADA, an International Federation, or a National Anti-Doping Organization*); or (2) another organization conducting *Testing* pursuant to the authority of and in accordance with the rules of the *Anti-Doping Organization* (for example, a National Federation that is a member of an International Federation).

## 8. Bibliography

1. *The World Anti-Doping Code International Standard for Laboratories* v8.0. World Anti-Doping Agency, Montreal, Canada (2015).  
<https://www.wada-ama.org/en/resources/laboratories/international-standard-for-laboratories-isl>
2. Powrie JK, Bassett EE, Rosen T *et al.* Detection of growth hormone abuse in sport. *Growth Horm IGF Res* **17**(3):220-6 (2007).
3. Guha N1, Erotokritou-Mulligan I, Bartlett C *et al.*. Biochemical markers of insulin-like growth factor-I misuse in athletes: the response of serum IGF-I, procollagen type III amino-terminal propeptide, and the GH-2000 score to the administration of rhIGF-I/rhIGF binding protein-3 complex.. *J Clin Endocrinol Metab.* **99**(6):2259-68 (2014). doi: 10.1210/jc.2013-3897. Epub 2014 Feb 25.
4. Holt RI, Sonksen PH. Growth hormone, IGF-I and insulin and their abuse in sport. *Br J Pharmacol* **154**(3):542-56 (2008).
5. Holt RI, Erotokritou-Mulligan I, McHugh C *et al.* The GH-2004 project: the response of IGF1 and type III pro-collagen to the administration of exogenous GH in non-Caucasian amateur athletes. *Eur J Endocrinol* **163**(1):45-54 (2010).
6. Erotokritou-Mulligan I, Bassett EE, Knies A, Sonksen PH, Holt RI. Validation of the growth hormone (GH)-dependent marker method of detecting GH abuse in sport through the use of independent data sets. *Growth Horm IGF Res* **17**(5):416-23 (2007).
7. Longobardi S, Keay N, Ehrnborg C *et al.* Growth hormone (GH) effects on bone and collagen turnover in healthy adults and its potential as a marker of GH abuse in sports: a double blind, placebo-controlled study. The GH-2000 Study Group. *J Clin Endocrinol Metab* **85**(4):1505-12 (2000).
8. Wallace JD, Cuneo RC, Lundberg PA *et al.* Responses of markers of bone and collagen turnover to exercise, growth hormone (GH) administration, and GH withdrawal in trained adult males. *J Clin Endocrinol Metab* **85**(1):124-33 (2000).
9. Wallace JD, Cuneo RC, Baxter R *et al.* Responses of the growth hormone (GH) and insulin-like growth factor axis to exercise, GH administration, and GH withdrawal in trained adult males: a potential test for GH abuse in sport. *J Clin Endocrinol Metab* **84**(10):3591-601 (1999).
10. Dall R, Longobardi S, Ehrnborg C *et al.* The effect of four weeks of supraphysiological growth hormone administration on the insulin-like growth factor axis in women and men. GH-2000 Study Group. *J Clin Endocrinol Metab* **85**(11):4193-200 (2000).
11. Nelson AE, Meinhardt U, Hansen JL *et al.* Pharmacodynamics of growth hormone abuse biomarkers and the influence of gender and testosterone: a randomized double-blind placebo-controlled study in young recreational athletes. *J Clin Endocrinol Metab* **93**(6):2213-22 (2008).
12. Erotokritou-Mulligan I, Bassett EE *et al.* Influence of ethnicity on IGF-I and procollagen III peptide (P-III-P) in elite athletes and its effect on the ability to detect GH abuse. *Clin Endocrinol (Oxf)* **70**(1):161-8 (2009).
13. Erotokritou-Mulligan I, Bassett E.E., Cowan DA *et al.* The use of growth hormone (GH)-dependent markers in the detection of GH abuse in sport: Physiological intra-individual variation of IGF-I, type 3 pro-collagen (P-III-P) and the GH-2000 detection score. *Clin Endocrinol (Oxf)* **72**(4):520-6 (2010).

14. Erotokritou-Mulligan I, Bassett EE, Bartlett C, Cowan D *et al.* The effect of sports injury on insulin-like growth factor-I and type 3 procollagen: implications for detection of growth hormone abuse in athletes. *J Clin Endocrinol Metab* **93**(7):2760-3 (2008).
15. Healy ML, Dall R, Gibney J *et al.* Toward the development of a test for growth hormone (GH) abuse: a study of extreme physiological ranges of GH-dependent markers in 813 elite athletes in the postcompetition setting. *J Clin Endocrinol Metab* **90**(2):641-9 (2005).
16. Nelson AE, Howe CJ, Nguyen TV *et al.* Influence of demographic factors and sport type on growth hormone-responsive markers in elite athletes. *J Clin Endocrinol Metab* **91**(11):4424-32 (2006).
17. Healy ML, Gibney J, Pentecost C, *et al.* Endocrine profiles in 693 elite athletes in the postcompetition setting. *Clin Endocrinol (Oxf)*. **81**(2):294-305 (2014). doi: 10.1111/cen.12445. Epub 2014 Apr 2.
18. Cox HD, Lopes F, Woldemariam GA *et al.* Interlaboratory Agreement of Insulin-like Growth Factor 1 Concentrations Measured by Mass Spectrometry. *Clinical Chemistry* **60**:541-548 (2014).
19. Bidlingmaier M, Friedrich N, Emeny RT *et al.* Reference intervals for insulin-like growth factor-1 (igf-i) from birth to senescence: results from a multicenter study using a new automated chemiluminescence IGF-I immunoassay conforming to recent international recommendations. *Clin Endocrinol Metab*. **99**: 1712-21 (2014).
20. Knudsen CS, Heickendorff L, Nexø E. Measurement of amino terminal propeptide of type III procollagen (PIIINP) employing the ADVIA Centaur platform. Validation, reference interval and comparison to UniQ RIA. *Clin Chem Lab Med* **52**(2): 237-241 (2014).
21. World Anti-Doping Program. Guidelines for Blood Sample Collection.  
[https://www.wada-ama.org/en/resources/search?f\[0\]=field\\_resource\\_collections%3A190](https://www.wada-ama.org/en/resources/search?f[0]=field_resource_collections%3A190)
22. Holt RI, Erotokritou-Mulligan I, Ridley SA *et al.* A determination of the pre-analytical storage conditions for insulin like growth factor-I and type III procollagen peptide. *Growth Horm IGF Res* **19**(1):43-50 (2009).
23. Holt RI, Böhning W, Guha N *et al.* The development of decision limits for the GH-2000 detection methodology using additional insulin-like growth factor-I and amino-terminal pro-peptide of type III collagen assays. *Drug Test Anal* **7**: 745-755 (2015).
24. WADA Technical Document TDDL: Decision Limits for the Confirmatory Quantification of Threshold Substances.  
[https://www.wada-ama.org/en/resources/search?f\[0\]=field\\_resource\\_collections%3A30](https://www.wada-ama.org/en/resources/search?f[0]=field_resource_collections%3A30)
25. Bassett EE, Erotokritou-Mulligan I. Statistical issues in implementing the marker method. *Growth Horm IGF Res* **19**(4):361-5 (2009).
26. ISO/IEC Guide 98-3: 2008. Evaluation of Measurement Data – Guide to the Expression of Uncertainty in Measurement (GUM) (2008).  
<http://www.bipm.org/en/publications/guides/gum.html>

## **Extended Bibliography describing the hGH Biomarkers Test**

- Abellan R, Ventura R, Palmi I *et al.* Immunoassays for the measurement of IGF-II, IGFBP-2 and -3, and ICTP as indirect biomarkers of recombinant human growth hormone misuse in sport. Values in selected population of athletes. *J Pharm Biomed Anal.* **48**(3):844-52 (2008).
- Abellan R, Ventura R, Pichini S *et al.* Effect of physical fitness and endurance exercise on indirect biomarkers of recombinant growth hormone misuse: insulin-like growth factor I and procollagen type III peptide. *Int J Sports Med.* **27**(12):976-83 (2006).
- Abellan R, Ventura R, Pichini S *et al.* Evaluation of immunoassays for the measurement of insulin-like growth factor-I and procollagen type III peptide, indirect biomarkers of recombinant human growth hormone misuse in sport. *Clin Chem Lab Med.* **43**(1):75-85 (2005).
- Armanini D, Faggian D, Scaroni C, Plebani M. Growth hormone and insulin-like growth factor I in a Sydney Olympic gold medalist. *Br J Sports Med.* **36**(2):148-9 (2002).
- Bassett EE, Erotokritou-Mulligan I. Statistical issues in implementing the marker method. *Growth Horm IGF Res.* **19**(4):361-5 (2009).
- Chung L, Baxter RC. Detection of growth hormone responsive proteins using SELDI-TOF mass spectrometry. *Growth Horm IGF Res.* **19**(4):383-7 (2009).
- Di Luigi L, Rigamonti AE, Agosti F *et al.* Combined evaluation of resting IGFI, N-terminal propeptide of type III procollagen and C-terminal cross-linked telopeptide of type I collagen levels might be useful for detecting inappropriate GH administration in female athletes. *Eur J Endocrinol.* **160**(5):753-8 (2009).
- Di Luigi L, Guidetti L. IGF-I, IGFBP-2, and -3: do they have a role in detecting rhGH abuse in trained men? *Med Sci Sports Exerc.* **34**(8):1270-8 (2002).
- Ding J, List EO, Okada S, Kopchick JJ. Perspective: proteomic approach to detect biomarkers of human growth hormone. *Growth Horm IGF Res.* **19**(4):399-407 (2009).
- Ehrnborg C, Lange KH, Dall R *et al*; GH-2000 Study Group. The growth hormone /insulin-like growth factor-I axis hormones and bone markers in elite athletes in response to a maximum exercise test. *J Clin Endocrinol Metab.* **88**(1):394-401 (2003).
- Erotokritou-Mulligan I, Eryl Bassett E, Cowan DA *et al.* The use of growth hormone (GH)-dependent markers in the detection of GH abuse in sport: Physiological intra-individual variation of IGF-I, type 3 pro-collagen (P-III-P) and the GH-2000 detection score. *Clin Endocrinol (Oxf).* **72**(4):520-6 (2010).
- Erotokritou-Mulligan I, Bassett EE, Cowan DA *et al*; GH-2004 group. Influence of ethnicity on IGF-I and procollagen III peptide (P-III-P) in elite athletes and its effect on the ability to detect GH abuse. *Clin Endocrinol (Oxf).* **70**(1):161-8 (2009).
- Erotokritou-Mulligan I, Bassett EE, Bartlett C *et al*; GH-2004 Group. The effect of sports injury on insulin-like growth factor-I and type 3 procollagen: implications for detection of growth hormone abuse in athletes. *J Clin Endocrinol Metab.* **93**(7):2760-3 (2008).
- Erotokritou-Mulligan I, Bassett EE, Knies A, Sonksen PH, Holt RI. Validation of the growth hormone (GH)dependent marker method of detecting GH abuse in sport through the use of independent data sets. *Growth Horm IGF Res.* **17**(5):416-23 (2007).
- Guha N, Erotokritou-Mulligan I, Burford C *et al.* Serum insulin-like growth factor-I and pro-collagen type III N-terminal peptide in adolescent elite athletes: implications for the detection of growth hormone abuse in sport. *J Clin Endocrinol Metab.* **95**(6):2969-76 (2010).
- Holt RI, Erotokritou-Mulligan I, McHugh C *et al.* The GH-2004 project: the response of IGF1 and type III pro-collagen to the administration of exogenous GH in non-Caucasian amateur athletes. *Eur J Endocrinol.* **163**(1):45-54 (2010).

- Holt RI, Bassett EE, Erotokritou-Mulligan I, McHugh C *et al*; GH-2004 group. Moving one step closer to catching the GH cheats: The GH-2004 experience. *Growth Horm IGF Res.* **19**(4):346-51 (2009).
- Holt RI, Erotokritou-Mulligan I, Sonksen PH. The history of doping and growth hormone abuse in sport. *Growth Horm IGF Res.* **19**(4):320-6 (2009).
- Kicman AT, Miell JP, Teale JD *et al*. Serum IGF-I and IGF binding proteins 2 and 3 as potential markers of doping with human GH. *Clin Endocrinol (Oxf).* **47**(1):43-50 (1997).
- Knies A, Ziegler E, Kratzsch J, Thieme D, Muller RK. Potential parameters for the detection of hGH doping. *Anal Bioanal Chem.* **376**(5):696-700 (2003).
- Longobardi S, Keay N, Ehrnborg C *et al*. Growth hormone (GH) effects on bone and collagen turnover in healthy adults and its potential as a marker of GH abuse in sports: a double blind, placebo-controlled study. The GH-2000 Study Group. *J Clin Endocrinol Metab.* **85**(4):1505-12 (2000).
- McHugh CM, Park RT, Sonksen PH, Holt RI. Challenges in detecting the abuse of growth hormone in sport. *Clin Chem.* **51**(9):1587-93 (2005).
- Minuto F, Barreca A, Melioli G. Indirect evidence of hormone abuse. Proof of doping? *J Endocrinol Invest.* **26**(9):919-23 (2003).
- Nelson AE, Ho KK. Demographic factors influencing the GH system: Implications for the detection of GH doping in sport. *Growth Horm IGF Res.* **19**(4):327-32 (2009).
- Nelson AE, Howe CJ, Nguyen TV *et al*. Influence of demographic factors and sport type on growth hormone-responsive markers in elite athletes. *J Clin Endocrinol Metab.* **91**(11):4424-32 (2006).
- Nguyen IV, Nelson AE, Howe CJ *et al*. Within subject variability and analytic imprecision of insulin like growth factor axis and collagen markers: implications for clinical diagnosis and doping tests. *Clin Chem.* **54**(8):1268-76 (2008).
- Pichini S, Ventura R, Palmi R *et al*. Effect of physical fitness and endurance exercise on indirect biomarkers of growth hormone and insulin misuse: Immunoassay-based measurement in urine samples. *J Pharm Biomed Anal.* **53**(4):1003-10 (2010).
- Powrie JK, Bassett EE, Rosen T *et al*; GH-2000 Project Study Group. Detection of growth hormone abuse in sport. *Growth Horm IGF Res.* **17**(3):220-6 (2007).
- Sartorio A, Jubeau M, Agosti F *et al*. A follow-up of GH dependent biomarkers during a 6-month period of the sporting season of male and female athletes. *J Endocrinol Invest.* **29**(3):237-43 (2006).
- Sartorio A, Agosti F, Marazzi N *et al*. Combined evaluation of resting IGF-I, N-terminal propeptide of type III procollagen (P-III-NP) and Cterminal cross-linked telopeptide of type I collagen (ICTP) levels might be useful for detecting inappropriate GH administration in athletes: a preliminary report. *Clin Endocrinol (Oxf).* **61**(4):487-93 (2004).
- Thevis M, Bredehöft M, Kohler M, Schänzer W. Mass spectrometry-based analysis of IGF-I and hGH. *Handb Exp Pharmacol.* **195**:201-7 (2010).
- Thevis M, Thomas A, Schänzer W. Doping control analysis of selected peptide hormones using LCMS(/MS). *Forensic Sci. Int.* **213**:35-41 (2011)
- Wallace JD, Cuneo RC, Baxter R *et al*. Responses of the growth hormone (GH) and insulin-like growth factor axis to exercise, GH administration, and GH withdrawal in trained adult males: a potential test for GH abuse in sport. *J Clin Endocrinol Metab.* **84**(10):3591-601 (1999).